

Anti-Atm Monoclonal Antibody For Wb, Ihc-P, Elisa - Q13315

Item Number: CM0002210



Introduction

ATM Rabbit monoclonal antibody validated for Western blot, IHC-P, and ELISA. Detects human and mouse ATM serine/threonine kinase, essential for DNA damage response and cell cycle control. Research use only.

[Learn More](#)

Parameter	Value
Product Name	ATM Rabbit mAb
Remarks/Alias	AT1, ATA, ATC, ATD, ATE, ATDC, TEL1, TELO1, ATM
Species	Human
GeneID (Human)	472
GeneID	472
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2619-2939 of human ATM (Q13315).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13315
Protein Weight	351kDa
Shipping	Ice bag