

Anti-Atm Rabbit Polyclonal Antibody For Wb, Elisa - Q13315

Item Number: CM0028845



Introduction

ATM Rabbit pAb validated for WB and ELISA. Detects human ATM, a 351 kDa serine/threonine kinase that senses DNA double-strand breaks and regulates cell cycle checkpoints. Ideal for research use.

[Learn More](#)

Parameter	Value
Product Name	ATM Rabbit pAb
Remarks/Alias	AT1, ATA, ATC, ATD, ATE, ATDC, TEL1, TELO1, ATM
Species	Human
GeneID (Human)	472
GeneID	472
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1900-2000 of human ATM (NP_000042.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13315
Protein Weight	351kDa
Shipping	Ice bag