

Anti-Atm Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13315

Item Number: CM0014127



Introduction

ATM Rabbit Polyclonal Antibody detects ATM (Q13315) in Human, Rat by WB, IF/ICC, ELISA. Serine/threonine kinase central to DNA damage response, checkpoint signaling, tumor suppression. Ideal for genomic stability research.

[Learn More](#)

Parameter	Value
Product Name	ATM Rabbit pAb
Remarks/Alias	AT1; ATA; ATC; ATD; ATE; ATDC; TEL1; TEL01; ATM
Species	Human
Gene ID (Human)	472
Gene ID	472
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 2320-2566 of ATM (NP_000042.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q13315
Protein Weight	351kDa
Shipping	Ice bag