

Anti-Atn1 Polyclonal Antibody For Wb, If/Icc, Elisa - P54259

Item Number: CM0021991

Introduction

High-quality rabbit polyclonal antibody against human ATN1 (atrophin-1), suitable for Western blot, immunofluorescence, and ELISA. Detects a 125 kDa protein involved in transcriptional regulation and neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	ATN1 Rabbit pAb
Remarks/Alias	B37, HRS, NOD, DRPLA, CHEDDA, D12S755E, ATN1
Species	Human
GeneID (human)	1822
GeneID	1822
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 900-1040 of human ATN1 (NP_001007027.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P54259
Protein Weight	125kDa
Shipping	Ice bag