

Anti-Atoh1 Monoclonal Antibody For Wb, Elisa - Q92858

Item Number: CM0003063

Introduction

High-specificity rabbit monoclonal antibody targeting human ATOH1 (MATH-1), validated for WB and ELISA. Ideal for neurogenesis and hearing loss research; detects a 38 kDa band. SWISS: Q92858.

[Learn More](#)

Parameter	Value
Product Name	ATOH1 Rabbit mAb
Remarks/Alias	ATH1; HATH1; DFNA89; MATH-1; bHLHa14; ATOH1
Species	Human
GeneID	474 (Human)
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ATOH1 (Q92858)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q92858
Protein Weight	38kDa
Shipping	Ice bag