

Anti-Atoh1 Polyclonal Antibody For Wb, Elisa - Q92858

Item Number: CM0028417

Introduction

Rabbit polyclonal antibody targeting human ATOH1 (bHLHa14). Validated for WB and ELISA, with cross-reactivity to human and mouse. Ideal for neurogenesis and transcriptional regulation studies. SWISS: Q92858.

[Learn More](#)

Parameter	Value
Product Name	ATOH1 Rabbit pAb
Remarks/Alias	ATH1; HATH1; DFNA89; MATH-1; bHLHa14; ATOH1
Species	Human
GeneID (Human)	474
GeneID	474
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 215-354 of human ATOH1 (NP_005163.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92858
Protein Weight	38kDa
Shipping	Ice bag