

Anti-Atoh8 Polyclonal Antibody For Wb And Elisa - Q96Sq7

Item Number: CM0022194

Introduction

Rabbit polyclonal antibody against human ATOH8 transcription factor. Validated for Western blot and ELISA applications. Suitable for research on endothelial cell biology, developmental processes, and gene regulation.

[Learn More](#)

Parameter	Value
Product Name	ATOH8 Rabbit pAb
Remarks/Alias	ATOH8; HATH6; bHLHa21; protein atonal homolog 8
Species	Human
GeneID	84913
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human ATOH8 (NP_116216.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96SQ7
Protein Weight	35kDa
Shipping	Ice bag