


Anti-Atox1 Rabbit Polyclonal Antibody For Wb, Elisa - O00244

Item Number: CM0029913



Introduction

Rabbit polyclonal antibody against human ATOX1 (copper transport protein ATX1). Validated for Western blot (WB) and ELISA. Suitable for copper homeostasis and antioxidant defense research.

[Learn More](#)

Parameter	Value
Product Name	ATOX1 Rabbit pAb
Remarks/Alias	ATX1; HAH1; ATOX1
Species	Human
GeneID (Human)	475
GeneID	475
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-68 of human ATOX1 (NP_004036.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O00244
Protein Weight	7kDa
Shipping	Ice bag