

Anti-Atp12A Polyclonal Antibody For Wb, Elisa - P54707

Item Number: CM0014133

Introduction

Rabbit anti-ATP12A polyclonal antibody validated for Western blot and ELISA in human samples. Detects the 116 kDa potassium-transporting ATPase alpha chain 2, involved in ion homeostasis and mucus acidification. For research use only.

[Learn More](#)

Parameter	Value
Product Name	ATP12A Rabbit pAb
Remarks/Alias	HK; ATP1AL1; H-K-ATPase; ATP12A
Species	Human
Gene ID (Human)	479
Gene ID	479
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 890-980 of human ATP12A (NP_001172014.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P54707
Protein Weight	116kDa
Shipping	Ice bag