


Anti-Atp1A4 Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q13733

Item Number: CM0007234



Introduction

Rabbit monoclonal antibody against human ATP1A4 (sperm-specific Na/K-ATPase alpha-4 subunit). Suitable for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. For research on sperm motility, capacitation, and fertility.

[Learn More](#)

Parameter	Value
Product Name	ATP1A4 Rabbit mAb
Remarks/Alias	ATP1A1; ATP1A2; ATP1A4
Species	Human
GeneID (human)	480
GeneID	480
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 401-530 of human ATP1A4 (NP_653300.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13733
Protein Weight	19KDa/114KDa
Shipping	Ice bag