


Anti-Atp1B1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P05026

Item Number: CM0005807



Introduction

Rabbit monoclonal antibody against human ATP1B1 (Sodium/potassium-transporting ATPase subunit beta-1). Validated for WB, IHC-P, ELISA. Cross-reacts with Human, Mouse, Rat. Suitable for studying ion transport, cell adhesion, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	ATP1B1 Rabbit mAb
Remarks/Alias	ATP1B; ATP1B1
Species	Human
GeneID (Human)	481
GeneID	481
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 204-303 of human ATP1B1 (P05026).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05026
Protein Weight	35kDa
Shipping	Ice bag