


## Anti-Atp1B1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P05026

Item Number: CM0025935



### Introduction

Rabbit polyclonal antibody targeting human ATP1B1 (sodium/potassium-transporting ATPase subunit beta-1), validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse and rat. Useful for studying ion transport, cell adhesion, and innate immunity.

[Learn More](#)

| Parameter        | Specification                                                                                                      |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | ATP1B1 Rabbit pAb                                                                                                  |
| Remarks/Alias    | ATP1B; ATP1B1                                                                                                      |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 481                                                                                                                |
| GeneID           | 481                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 63-303 of human ATP1B1 (NP_001668.1) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | P05026                                                                                                             |
| Protein Weight   | 35kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |