

Anti-Atp1B2 Polyclonal Antibody For Ihc-P, Elisa - P14415

Item Number: CM0019567

Introduction

Rabbit polyclonal antibody targeting human ATP1B2 (Sodium/potassium-transporting ATPase subunit beta-2). Validated for IHC-P and ELISA, with mouse and rat cross-reactivity. Ideal for cell adhesion and ion transport research.

[Learn More](#)

Parameter	Value
Product Name	ATP1B2 Rabbit pAb
Remarks/Alias	AMOG; ATP1B2
Species	Human
GeneID (Human)	482
GeneID	482
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 68-290 of human ATP1B2 (NP_001669.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P14415
Protein Weight	33kDa
Shipping	Ice bag