

Anti-Atp1B2 Rabbit Polyclonal Antibody For Wb - P14415

Item Number: CM0028506

Introduction

High-quality rabbit polyclonal antibody against ATP1B2 (AMOG), validated for WB, IF/ICC, ELISA. Reacts with human, mouse, and rat samples. Ideal for IVD research on sodium/potassium ATPase function and cell adhesion mechanisms.

[Learn More](#)

Parameter	Value
Product Name	ATP1B2 Rabbit pAb
Remarks/Alias	AMOG; ATP1B2
Species	Human
Gene ID (Human)	482
Gene ID	482
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 68-290 of human ATP1B2 (NP_001669.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14415
Protein Weight	33kDa
Shipping	Ice bag