


Anti-Atp1B3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P54709

Item Number: CM0003217



Introduction

High-quality ATP1B3 rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Recognizes human, mouse, and rat ATP1B3 (CD298/ATPB-3) with observed bands at 22 and 32 kDa. Ideal for ion transport and membrane protein research.

[Learn More](#)

Parameter	Details
Product Name	ATP1B3 Rabbit mAb
Remarks/Alias	CD298; ATPB-3; ATP1B3
Species	Human
GeneID (Human)	483
GeneID	483
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 180-279 of human ATP1B3 (P54709)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54709
Protein Weight	22kDa/32kDa
Shipping	Ice bag