

Anti-Atp23 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y6H3

Item Number: CM0017831

Introduction

Rabbit polyclonal antibody against ATP23 (KUB3/XRCC6BP1) for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for mitochondrial protease and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	ATP23 Rabbit pAb
Remarks/Alias	KUB3; XRCC6BP1; ATP23
Species	Human
Gene ID (Human)	91419
Gene ID	91419
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-246 of human ATP23 (NP_150592.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6H3
Protein Weight	28kDa
Shipping	Ice bag