

Anti-Atp2B1 Monoclonal Antibody For Wb, If/Icc, Elisa - P20020

Item Number: CM0006344



Introduction

Anti-ATP2B1 rabbit monoclonal antibody, validated in WB, IF/ICC, ELISA, cross-reacts with human, mouse, rat. For studying calcium homeostasis, blood pressure regulation, bone mineralization, and insulin sensitivity.

[Learn More](#)

Parameter	Value
Product Name	ATP2B1 Rabbit mAb
Remarks/Alias	MRD66; PMCA1; PMCA1kb; ATP2B1
Species	Human
GeneID	490
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20020
Protein Weight	135 kDa
Shipping	Ice bag