

Anti-Atp2B1 Polyclonal Antibody For Wb, If-P, Elisa - P20020

Item Number: CM0021516



Introduction

The ATP2B1 Rabbit Polyclonal Antibody (pAb) detects ATP2B1/PMCA1 in human, mouse, rat samples. Validated for WB, IF-P, and ELISA; ideal for studying calcium transport, blood pressure regulation, and bone metabolism.

[Learn More](#)

Parameter	Value
Product Name	ATP2B1 Rabbit pAb
Remarks/Alias	MRD66; PMCA1; PMCA1kb; ATP2B1
Species	Human
Gene ID (Human)	490
Gene ID	490
Immunogen	Synthetic peptide corresponding to amino acids 50-150 of human ATP2B1 (NP_001673.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20020
Protein Weight	135kDa
Shipping	Ice bag