

Anti-Atp2B2 Rabbit Polyclonal Antibody For Wb, Elisa - Q01814

Item Number: CM0028421

Introduction

High-quality rabbit polyclonal antibody targeting ATP2B2 (PMCA2), suitable for WB and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying calcium homeostasis and auditory function.

[Learn More](#)

Parameter	Value
Product Name	ATP2B2 Rabbit pAb
Remarks/Alias	PMCA2; DFNA82; PMCA2a; PMCA2i; ATP2B2
Species	Human
GeneID (Human)	491
GeneID	491
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human ATP2B2 (NP_001674.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01814
Protein Weight	137kDa
Shipping	Ice bag