

Anti-Atp4B (Hydrogen Potassium Atpase Beta) Rabbit Monoclonal Antibody For Ihc-P, Elisa - P51164

Item Number: CM0003296

Introduction

Rabbit monoclonal antibody targeting human ATP4B/Potassium-transporting ATPase subunit beta. Validated for IHC-P and ELISA; cross-reacts with mouse and rat. Supports gastric proton pump and ion transport studies.

[Learn More](#)

Parameter	Value
Product Name	ATP4B Rabbit mAb
Remarks/Alias	Hydrogen Potassium ATPase Beta; ATP4B
Species	Human
GeneID (Human)	496
GeneID	496
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 192-291 of human ATP4B (P51164)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51164
Protein Weight	33kDa
Shipping	Ice bag