


Anti-Atp5A1 Polyclonal Antibody For Wb, If/Icc, Elisa - P25705

Item Number: CM0029425



Introduction

Validated in WB, IF/ICC, ELISA, this rabbit polyclonal antibody targets human ATP5A1 (ATP synthase F1 subunit alpha), cross-reacting with mouse and rat. Ideal for mitochondrial research and energy metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	ATP5A1 Rabbit pAb
Remarks/Alias	OMR; ORM; ATPM; MOM2; ATP5A; hATP1; ATP5A1; MC5DN4; ATP5AL2; COXPD22; HEL-S-123m
Species	Human
Gene ID (Human)	498
Gene ID	498
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human ATP5A1 (NP_004037.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25705
Protein Weight	60kDa
Shipping	Ice bag