


Anti-Atp5A1 Rabbit Monoclonal Antibody For Wb, Ip, If, Ihc-P, Elisa - P25705

Item Number: CM0002190



Introduction

Rabbit monoclonal antibody against ATP5A1, validated for WB, IP, IF/ICC, IHC-P, and ELISA. Detects human, mouse, and rat ATP synthase subunit alpha. Ideal for mitochondrial function studies.

[Learn More](#)

Parameter	Value
Product Name	ATP5A1 Rabbit mAb
Remarks/Alias	OMR; ORM; ATPM; MOM2; ATP5A; hATP1; ATP5A1; MC5DN4; ATP5AL2; COXPD22; HEL-S-123m
Species	Human
Gene ID (Human)	498
Gene ID	498
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human ATP5A1 (P25705).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25705
Protein Weight	60kDa
Shipping	Ice bag