


Anti-Atp5C1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P36542

Item Number: CM0013920



Introduction

ATP5C1 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat ATP synthase gamma subunit (P36542), a central stalk component of mitochondrial Complex V, suitable for metabolism and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	ATP5C1 Rabbit pAb
Remarks/Alias	ATP5C; ATP5C1; ATP5CL1
Species	Human
Gene ID	509
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-298 of human ATP5C1 (NP_001001973.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36542
Protein Weight	33kDa
Shipping	Ice bag