

Anti-Atp5D Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P30049

Item Number: CM0026214



Introduction

Rabbit polyclonal antibody targeting human ATP5D (ATP synthase F1 subunit delta, mitochondrial). Validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for mitochondrial research and ATP synthase studies.

[Learn More](#)

Parameter	Description
Product Name	ATP5D Rabbit pAb
Remarks/Alias	ATP5D; MC5DN5
Species	Human
GeneID (human)	513
GeneID	513
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-168 of human ATP5D (NP_001678.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30049
Protein Weight	17kDa
Shipping	Ice bag