

Anti-Atp5F1 Polyclonal Antibody For Wb, If/Icc, Elisa - P24539

Item Number: CM0025804



Introduction

High-quality rabbit polyclonal antibody targeting ATP5F1, a subunit of mitochondrial ATP synthase. Validated for WB, IF/ICC, ELISA; reacts with human and mouse. Ideal for mitochondrial bioenergetics research.

[Learn More](#)

Parameter	Value
Product Name	ATP5F1 Rabbit pAb
Remarks/Alias	PIG47; ATP5F1
Species	Human
GeneID (Human)	515
GeneID	515
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-245 of human ATP5F1 (NP_001679.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P24539
Protein Weight	29kDa
Shipping	Ice bag