

Anti-Atp5G2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q06055

Item Number: CM0020479

Introduction

Rabbit polyclonal antibody against human ATP5G2, validated for WB, IF/ICC, and ELISA. Reacts with human and mouse samples. Ideal for mitochondrial ATP synthase research and proton transport studies.

[Learn More](#)

Parameter	Value
Product Name	ATP5G2 Rabbit pAb
Remarks/Alias	ATP5A, ATP5G2
Species	Human
Gene ID (Human)	517
Gene ID	517
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human ATP5G2 (NP_001002031.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q06055
Protein Weight	15kDa
Shipping	Ice bag