

Anti-Atp5H Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O75947

Item Number: CM0028904



Introduction

High-quality rabbit polyclonal antibody against ATP5H (ATP synthase peripheral stalk subunit d, O75947). Validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for mitochondrial metabolism research.

[Learn More](#)

Parameter	Value
Product Name	ATP5H Rabbit pAb
Remarks/Alias	ATPQ; APT5H; ATP5H
Species	Human
GeneID (Human)	10476
GeneID	10476
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-161 of human ATP5H (NP_006347.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75947
Protein Weight	18kDa
Shipping	Ice bag