

Anti-Atp5I Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P56385

Item Number: CM0020709



Introduction

Rabbit polyclonal antibody targeting human ATP5I (ATP synthase subunit e). Validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for mitochondrial and ATP synthase studies.

[Learn More](#)

Parameter	Value
Product Name	ATP5I Rabbit pAb
Remarks/Alias	ATP5I; ATP5K
Species	Human
GeneID (Human)	521
GeneID	521
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56385
Protein Weight	8kDa
Shipping	Ice bag