

Anti-Atp5j2 Polyclonal Antibody For Wb And Elisa - P56134

Item Number: CM0019905



Introduction

Rabbit polyclonal antibody targeting human ATP5J2, suitable for WB and ELISA applications, cross-reacts with mouse and rat. Ideal for mitochondrial ATP synthase complex studies.

[Learn More](#)

Parameter	Value
Product Name	ATP5J2 Rabbit Polyclonal Antibody
Remarks/Alias	ATP5J2, ATP5JL
Species	Human
Gene ID	9551
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-94 of human ATP5J2 (NP_004880.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P56134
Protein Weight	11kDa
Shipping	Ice bag