

## Anti-Atp5L Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O75964

Item Number: CM0026159



### Introduction

Rabbit polyclonal antibody against human ATP5L (ATP synthase membrane subunit g), suitable for WB, IF/ICC, ELISA. Reacts with human, mouse, and rat. Ideal for mitochondrial ATP synthase research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | ATP5L Rabbit pAb                                                                                                  |
| Remarks/Alias    | ATP5L; ATP5JG                                                                                                     |
| Species          | Human                                                                                                             |
| Gene ID (Human)  | 10632                                                                                                             |
| Gene ID          | 10632                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-103 of human ATP5L (NP_006467.4). |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, IF/ICC, ELISA                                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | O75964                                                                                                            |
| Protein Weight   | 11kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |