

Anti-Atp6 Rabbit Monoclonal Antibody For Wb, Elisa - P00848

Item Number: CM0006710

Introduction

Rabbit monoclonal antibody against ATP6 (ATP synthase subunit a) for WB and ELISA. Detects human ATP6, a key proton channel subunit of mitochondrial ATP synthase. Ideal for metabolism and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	ATP6 Rabbit mAb
Remarks/Alias	ATP6
Species	Mouse
GeneID (Human)	4508
GeneID	17705
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P00848
Protein Weight	25kDa
Shipping	Ice bag