


Anti-Atp6Ap1 Polyclonal Antibody For Wb, Elisa - Q15904

Item Number: CM0017529



Introduction

Rabbit polyclonal antibody against ATP6AP1 (V-ATPase subunit S1), validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Ideal for studying proton pump function, membrane trafficking, and intracellular iron homeostasis.

[Learn More](#)

Parameter	Value
Product Name	ATP6AP1 Rabbit pAb
Remarks/Alias	16A; CF2; Ac45; XAP3; XAP-3; ATP6S1; VATPS1; ATP6IP1; ATP6AP1
Species	Human
GeneID (Human)	537
GeneID	537
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-390 of human ATP6AP1 (NP_001174.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15904
Protein Weight	52kDa
Shipping	Ice bag