


Anti-Atp6Ap1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q15904

Item Number: CM0006821



Introduction

Recombinant rabbit monoclonal antibody targeting human ATP6AP1 (V-ATPase subunit S1, ~52 kDa). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse. Ideal for lysosomal acidification and membrane trafficking research.

[Learn More](#)

Parameter	Value
Product Name	ATP6AP1 Rabbit mAb
Remarks/Alias	16A; CF2; Ac45; XAP3; XAP-3; ATP6S1; VATPS1; ATP6IP1; ATP6AP1
Species	Human
GeneID (Human)	537
GeneID	537
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15904
Protein Weight	52kDa
Shipping	Ice bag