


Anti-Atp6Ap2 Rabbit Polyclonal Antibody For Wb, Elisa - 075787

Item Number: CM0026735



Introduction

Anti-ATP6AP2 Rabbit Polyclonal Antibody validated for WB and ELISA. Recognizes human, mouse, rat ATP6AP2 protein (renin receptor, 39 kDa). Suitable for IVD research in renin-angiotensin system and lysosomal biology.

[Learn More](#)

Parameter	Value
Product Name	ATP6AP2 Rabbit pAb
Remarks/Alias	PRR; M8-9; MRXE; RENR; XMRE; XPDS; CDG2R; HT028; MRXSH; ELDF10; ATP6IP2; MSTP009; APT6M8-9; ATP6M8-9; ATP6AP2
Species	Human
GeneID (Human)	10159
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 251-350 of human ATP6AP2 (NP_005756.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075787
Protein Weight	39kDa
Shipping	Ice bag