


Anti-Atp6V0A1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q93050

Item Number: CM0015621



Introduction

ATP6V0A1 Rabbit Polyclonal Antibody for WB, IF/ICC, ELISA; recognizes human, mouse, rat. Targets V-ATPase subunit a1 for organelle acidification, mTORC1 signaling, and neuronal development studies.

[Learn More](#)

Parameter	Value
Product Name	ATP6V0A1 Rabbit pAb
Remarks/Alias	a1; Stv1; VPP1; Vph1; ATP6N1; DEE104; NEDEBA; ATP6N1A; ATP6V0A1
Species	Human
GeneID (human)	535
GeneID	535
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-300 of human ATP6V0A1 (NP_001123493.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q93050
Protein Weight	96kDa
Shipping	Ice bag