


Anti-Atp6V0C Rabbit Monoclonal Antibody For Wb And Elisa - P27449

Item Number: CM0004716



Introduction

High-quality rabbit monoclonal antibody targeting human ATP6V0C (V-type proton ATPase 16 kDa proteolipid subunit c), validated for Western blot and ELISA applications. Ideal for studies on vacuolar acidification and lysosomal function.

[Learn More](#)

Parameter	Value
Product Name	ATP6V0C Rabbit mAb
Remarks/Alias	ATP6V0C; ATP6C; ATP6L; ATPL; VATL; VPPC; Vma3; ATPase H ⁺ transporting V0 subunit c
Species	Human
Gene ID (Human)	527
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ATP6V0C. (NP_001685.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P27449
Protein Weight	15kDa
Shipping	Ice bag