


Anti-Atp6V0C Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P27449

Item Number: CM0023098



Introduction

Anti-ATP6V0C Rabbit Polyclonal Antibody, validated for Western blot (WB), Immunofluorescence/Immunocytochemistry (IF/ICC), and ELISA. Cross-reacts with Human, Mouse, and Rat. Essential for studying vacuolar ATPase proton channel function and lysosomal acidification.

[Learn More](#)

Parameter	Value
Product Name	ATP6V0C Rabbit pAb
Remarks/Alias	ATPL; VATL; VPPC; Vma3; ATP6C; ATP6L; ATP6V0C
Species	Human
GeneID (Human)	527
GeneID	527
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ATP6V0C (NP_001185498.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27449
Protein Weight	16kDa
Shipping	Ice bag