


Anti-Atp6V0D1 Rabbit Polyclonal Antibody For Wb, Elisa - P61421

Item Number: CM0031029



Introduction

Rabbit polyclonal antibody targeting human ATP6V0D1 (V-ATPase subunit d1), validated for WB and ELISA. Reacts with human, mouse, and rat. Ideal for studying lysosomal acidification and mTORC1 signaling.

[Learn More](#)

Parameter	Value
Product Name	ATP6V0D1 Rabbit pAb
Remarks/Alias	P39; VATX; VMA6; ATP6D; ATP6DV; VPATPD; ATP6V0D1
Species	Human
GeneID (human)	9114
GeneID	9114
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-351 of human ATP6V0D1 (NP_004682.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61421
Protein Weight	40kDa
Shipping	Ice bag