


Anti-Atp6V1A Polyclonal Antibody For Wb And Elisa - P38606

Item Number: CM0013720



Introduction

Rabbit polyclonal antibody targeting human ATP6V1A (V-ATPase catalytic subunit A). Validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for V-ATPase, ion transport, and autophagy research.

[Learn More](#)

Parameter	Value
Product Name	ATP6V1A Rabbit pAb
Remarks/Alias	HO68; VA68; VPP2; Vma1; DEE93; ARCL2D; ATP6A1; IECEE3; ATP6V1A1; ATP6V1A
Species	Human
GeneID (Human)	523
GeneID	523
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 500-600 of human ATP6V1A (NP_001681.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P38606
Protein Weight	68kDa
Shipping	Ice bag