

Anti-Atp6V1B1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P15313

Item Number: CM0030681

Introduction

High-quality rabbit polyclonal antibody targeting human ATP6V1B1 (V-ATPase subunit B1). Validated for WB, IF-P, ELISA with cross-reactivity to mouse and rat. Ideal for renal acidification and deafness research.

[Learn More](#)

Parameter	Value
Product Name	ATP6V1B1 Rabbit pAb
Remarks/Alias	VATB; VMA2; VPP3; DRTA2; RTA1B; ATP6B1; ATP6V1B1
Species	Human
Gene ID (Human)	525
Gene ID	525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human ATP6V1B1 (NP_001683.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15313
Protein Weight	57kDa
Shipping	Ice bag