


Anti-Atp6V1B2 Rabbit Polyclonal Antibody For Wb, Elisa - P21281

Item Number: CM0026824



Introduction

Rabbit polyclonal antibody against ATP6V1B2, suitable for Western blot and ELISA. Recognizes human, mouse, and rat V-ATPase subunit B brain isoform, a 57 kDa protein essential for proton transport and lysosomal acidification.

[Learn More](#)

Parameter	Value
Product Name	ATP6V1B2 Rabbit pAb
Remarks/Alias	DOOD; HO57; VATB; VPP3; Vma2; ZLS2; ATP6B2; ATP6B1B2; ATP6V1B2
Species	Human
GeneID (Human)	526
GeneID	526
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 262-511 of human ATP6V1B2 (NP_001684.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21281
Protein Weight	57kDa
Shipping	Ice bag