


Anti-Atp6V1D Polyclonal Antibody For Wb, Elisa - Q9Y5K8

Item Number: CM0015206



Introduction

Rabbit polyclonal antibody targeting human ATP6V1D (V-ATPase subunit D), suitable for WB and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying vacuolar acidification, mTORC1 signaling, and cilium biogenesis.

[Learn More](#)

Parameter	Value
Product Name	ATP6V1D Rabbit pAb
Remarks/Alias	VATD; VMA8; ATP6M; ATP6V1D
Species	Human
Gene ID (Human)	51382
Gene ID	51382
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-247 of human ATP6V1D (NP_057078.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5K8
Protein Weight	28kDa
Shipping	Ice bag