


Anti-Atp6V1E1 Polyclonal Antibody For Wb, If/Icc, Elisa - P36543

Item Number: CM0026234



Introduction

High-quality rabbit polyclonal antibody targeting ATP6V1E1, suitable for Western blot, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Essential for V-ATPase complex and ion transport research. - P36543

[Learn More](#)

Parameter	Value
Product Name	ATP6V1E1 Rabbit pAb
Remarks/Alias	P31; Vma4; ATP6E; ARCL2C; ATP6E2; ATP6V1E; ATP6V1E1
Species	Human
Gene ID (human)	529
Gene ID	529
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 77-226 of human ATP6V1E1 (NP_001687.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36543
Protein Weight	26kDa
Shipping	Ice bag