

Anti-Atp6V1G2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O95670

Item Number: CM0023654

Introduction

CamelBio's anti-ATP6V1G2 rabbit polyclonal antibody detects V-ATPase subunit G2 (14 kDa) in human, mouse, rat samples via WB, IHC-P, ELISA for vacuolar proton pump research.

[Learn More](#)

Parameter	Specification
Product Name	ATP6V1G2 Rabbit pAb
Remarks/Alias	NG38; ATP6G; VMA10; ATP6G2; ATP6V1G2
Species	Human
GeneID (Human)	534
GeneID	534
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-100 of human ATP6V1G2 (NP_569730.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95670
Protein Weight	14kDa
Shipping	Ice bag