

## Anti-Atp6V1G3 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96Lb4

Item Number: CM0016997

### Introduction

CamelBio offers a rabbit polyclonal antibody against ATP6V1G3, targeting the V-ATPase subunit G 3. Validated for WB, IHC-P, and ELISA with cross-reactivity to mouse. Ideal for renal acidification and ion transport research.

[Learn More](#)

| Parameter        | Value                                                                                                                                     |
|------------------|-------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | ATP6V1G3 Rabbit pAb                                                                                                                       |
| Remarks/Alias    | Vma10; ATP6G3; ATP6V1G3                                                                                                                   |
| Species          | Human                                                                                                                                     |
| GeneID (Human)   | 127124                                                                                                                                    |
| GeneID           | 127124                                                                                                                                    |
| Immunogen        | Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-118 of human ATP6V1G3 (NP_573569.1). |
| Source           | Rabbit                                                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                                                     |
| Application      | WB, IHC-P, ELISA                                                                                                                          |
| Cross Reactivity | Mouse                                                                                                                                     |
| SWISS            | Q96LB4                                                                                                                                    |
| Protein Weight   | 14kDa                                                                                                                                     |
| Shipping         | Ice bag                                                                                                                                   |