

Anti-Atp7A Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q04656

Item Number: CM0008382



Introduction

High-quality ATP7A rabbit monoclonal antibody validated for human sample applications including WB, IF/ICC, and ELISA. Ideal for investigating copper homeostasis, Menkes disease, and protein trafficking.

[Learn More](#)

Parameter	Value
Product Name	ATP7A Rabbit mAb
Remarks/Alias	MK; MNK; DSMAX; SMAX3
Species	Human
GeneID (Human)	538
GeneID	538
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-79 of human ATP7A (NP_000043.4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q04656
Protein Weight	11-172 kDa
Shipping	Ice bag