

Anti-Atp7A Rabbit Polyclonal Antibody For Wb And Elisa - Q04656

Item Number: CM0027360



Introduction

High specificity anti-ATP7A rabbit polyclonal antibody validated for WB and ELISA. Detects human ATP7A with cross-reactivity to rat. Ideal for copper homeostasis, Menkes disease, and neuronal function research.

[Learn More](#)

Parameter	Value
Product Name	ATP7A Rabbit pAb
Remarks/Alias	MK; MNK; DSMAX; SMAX3; ATP7A
Species	Human
Gene ID (Human)	538
Gene ID	538
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-79 of human ATP7A (NP_000043.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q04656
Protein Weight	11-172 kDa
Shipping	Ice bag