

Anti-Atp7B Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P35670

Item Number: CM0029077



Introduction

Rabbit polyclonal antibody targeting human copper-transporting ATPase 2 (ATP7B). Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for Wilson disease and copper homeostasis studies.

[Learn More](#)

Parameter	Specification
Product Name	ATP7B Rabbit pAb
Remarks/Alias	WD; PWD; WC1; WND; ATP7B
Species	Human
GeneID (Human)	540
GeneID	540
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35670
Protein Weight	157 kDa
Shipping	Ice bag