

Anti-Atp8 Monoclonal Antibody For Wb, Ip, Ihc-P, Elisa - P03928

Item Number: CM0005264



Introduction

High-quality ATP8 rabbit monoclonal antibody (P03928) validated for Western blot, IP, IHC-P, and ELISA. Reacts with human and rat ATP synthase F(0) complex subunit 8, with a molecular weight of 8 kDa.

[Learn More](#)

Parameter	Value
Product Name	ATP8 Rabbit mAb
Remarks/Alias	ATPase8; MTATP8; ATP8
Species	Human
GeneID (human)	4509
GeneID	4509
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 1-68 of human ATP8 (YP_003024030.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P03928
Protein Weight	8kDa
Shipping	Ice bag