

## Anti-Atp8B3 Rabbit Polyclonal Antibody For Wb And Elisa - O60423

Item Number: CM0020605

### Introduction

Rabbit polyclonal antibody against ATP8B3, a phospholipid-transporting ATPase flippase critical for sperm function. Validated for WB and ELISA with human samples. Protein weight 147kDa.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | ATP8B3 Rabbit pAb                                                                                                   |
| Remarks/Alias    | ATPIK; ATP8B3                                                                                                       |
| Species          | Human                                                                                                               |
| Gene ID (Human)  | 148229                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 600-900 of human ATP8B3 (NP_620168.1) |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, ELISA                                                                                                           |
| Cross Reactivity | Human                                                                                                               |
| SWISS            | O60423                                                                                                              |
| Protein Weight   | 147kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                             |